


# CRYOcheck™ CorPac™

CRYOcheck CorPac consists of 30 individual vials containing human plasma with distinct PT, APTT and fibrinogen profiles spanning the normal and pathological range. Plasmas are buffered

and rapidly frozen and are available in sets of 30 x 1.5 mL aliquots. CRYOcheck CorPac is intended *For Research Use Only*.

## FEATURES:

- Donors are screened negative for lupus anticoagulants
- Source plasmas are screened negative for all FDA-required tests
- Packaging and vials are compact and color-coded for easy storage and location in freezers
- Convenient frozen format – ready to use within minutes of thawing
- Expiry dating of 2 years from the date of manufacture when stored at -40 to -80 °C

DESCRIPTION	CATALOG #	FORMAT
 CRYOcheck CorPac	CCCP-15	30 x 1.5 mL



PrecisionBioLogic

[www.precisionbiologic.com](http://www.precisionbiologic.com)  
1.800.267.2796

Frozen diagnostic products. One simple idea, many cool benefits.