

CRYOcheck™ dRVVT Reagents



CRYOcheck LA Check™ (dRVVT Screening Reagent) and CRYOcheck LA Sure™ (dRVVT Confirmatory Reagent) are frozen dilute Russell's Viper Venom Time (dRVVT) reagents. These products can be integrated into a lupus anticoagulant (LA) test panel to provide users with enhanced convenience and increased confidence in testing results.

FEATURES:

- No reconstitution errors
- Easily adapted to manual, semi-automated, and automated instrumentation (protocols for most instrumentation available upon request)
- Unused reagent may be refrozen once
- Packaging and vials are compact and color-coded for easy storage and location in freezers
- Convenient frozen format – ready to use within minutes of thawing
- 48-hour stability once thawed when maintained at 2 to 8 °C
- Expiry dating of 3 years from the date of manufacture when stored at -40 to -80 °C
- Screening and confirmatory reagents are sold separately to eliminate excess stock of either product in your lab

CRYOcheck LA Check and CRYOcheck LA Sure are part of Precision BioLogic's suite of LA products. They can be used in conjunction with our controls: Lupus Positive Control, Weak Lupus Positive Control, Lupus Negative Control, and Pooled Normal Plasma. For a complete LA testing panel, Precision BioLogic offers Platelet Lysate for use in the platelet neutralization procedure and Hex LA,™ an integrated screen and confirm silica-based APTT assay, to aid in the detection of LA by the application of hexagonal phase phospholipids.

CRYOcheck LA Check is a sensitive reagent that detects weak LA samples. Reagent contains Russell's viper venom, phospholipids and anti-heparin agents. CRYOcheck LA Sure has a high phospholipid concentration that neutralizes LA and therefore demonstrates the correction needed to report positive LA results with confidence. The result is top quality products which can be used as part of a complete LA workup.

	DESCRIPTION	CATALOG #	FORMAT
	LA Check	CHK-10 CHK80-10	25 x 1.0 mL 80 x 1.0 mL
	LA Sure	SUR-10	25 x 1.0 mL



PrecisionBioLogic

www.precisionbiologic.com
1.800.267.2796

Frozen diagnostic products. One simple idea, many cool benefits.