

CRYOcheck™ Factor Deficient Plasmas

CRYOcheck Factor Deficient Plasmas are intended for use as deficient substrates in clot-based factor assays using the one-stage prothrombin time (PT) and activated partial thromboplastin time (APTT), and provide users with enhanced convenience and increased confidence in testing results. CRYOcheck Factor Deficient Plasmas










are frozen, platelet-poor plasmas which are immunodepleted and assayed at less than 1% for the specific factor by functional and antigenic methods*. The plasma is buffered and frozen, resulting in top quality products that are consistent vial to vial.

FEATURES:

- No reconstitution errors
- No bovine additives
- Immunodepleted factor deficient plasmas are screened negative for the presence of factor-specific inhibitors
- Packaging and vials are compact and color-coded for easy storage and location in freezers
- Convenient frozen format – ready to use within minutes of thawing
- 8-hour stability once thawed when maintained at 2 to 8 °C
- Expiry dating of 3 years from the date of manufacture when stored at -40 to -80 °C
- Quality Control Certificate provided for each lot number
- Plasma components are screened negative for all FDA-required tests
- Comprehensive selection: Factors II through XII and Prekallikrein

CRYOcheck Factor Deficient Plasmas may be used in conjunction with the members of our CRYOcheck Gold Standard Family of Assayed Plasmas, including Normal Reference Plasma, Reference Control Normal and Abnormal Reference Controls.

*CRYOcheck Prekallikrein Deficient Plasma is prepared from a prekallikrein-deficient donor.

DESCRIPTION	CATALOG #	FORMAT
 Factor II Deficient Plasma	FDP02-10 FDP02-15	25 x 1.0 mL 25 x 1.5 mL
 Factor V Deficient Plasma	FDP05-10 FDP05-15	25 x 1.0 mL 25 x 1.5 mL
 Factor VII Deficient Plasma	FDP07-10 FDP07-15	25 x 1.0 mL 25 x 1.5 mL
 Factor VIII Deficient Plasma	FDP08-10 FDP08-15	25 x 1.0 mL 25 x 1.5 mL
 Factor IX Deficient Plasma	FDP09-10 FDP09-15	25 x 1.0 mL 25 x 1.5 mL
 Factor X Deficient Plasma	FDP10-10 FDP10-15	25 x 1.0 mL 25 x 1.5 mL
 Factor XI Deficient Plasma	FDP11-10 FDP11-15	25 x 1.0 mL 25 x 1.5 mL
 Factor XII Deficient Plasma	FDP12-10 FDP12-15	25 x 1.0 mL 25 x 1.5 mL
 Prekallikrein Deficient Plasma	FDPK-10	10 x 1.0 mL



PrecisionBioLogic

www.precisionbiologic.com
1.800.267.2796

Frozen diagnostic products. One simple idea, many cool benefits.