

## CRYOcheck™ Gold Standard Family of Assayed Plasmas

Our Gold Standard Family of Assayed Plasmas is comprised of CRYOcheck Normal Reference Plasma (calibration plasma), CRYOcheck Reference Control Normal (normal control), CRYOcheck Abnormal 1 Reference Control (borderline pathological control), and CRYOcheck Abnormal 2 Reference Control (pathological control). These products are intended for use as assayed calibration and quality control plasmas for your special coagulation needs, and provide users with enhanced convenience and increased confidence in testing results.





Members of our Gold Standard Family of Assayed Plasmas are assayed for over 20 hemostatic parameters against international

reference standards (where available). Reference values are assigned by independent, internationally recognized reference laboratories. Carefully screened plasma from male and female donors between 18 and 66 years of age are used to manufacture all products in our Gold Standard Family of Assayed Plasmas. The plasma is buffered and frozen, resulting in top quality products that, together, provide assurance in the integrity of test results.

### FEATURES:

- No reconstitution errors
- 8-hour stability once thawed when maintained at 2 to 8 °C
- Packaging and vials are compact and color-coded for easy storage and location in freezers
- Convenient frozen format – ready to use within minutes of thawing
- Expiry dating of 3 years from the date of manufacture when stored at -40 to -80 °C
- Plasma components are screened negative for all FDA-required tests

Members of our Gold Standard Family of Assayed Plasmas may be used in conjunction with CRYOcheck Factor Deficient Plasmas.

DESCRIPTION	CATALOG #	FORMAT
 Normal Reference Plasma	CCNRP-05	25 x 0.5 mL
	CCNRP-10	25 x 1.0 mL
 Reference Control Normal	RCN-05	25 x 0.5 mL
	RCN-10	25 x 1.0 mL
 Abnormal 1 Reference Control	ARP1-05	25 x 0.5 mL
	ARP1-10	25 x 1.0 mL
 Abnormal 2 Reference Control	ARP2-10	25 x 1.0 mL



PrecisionBioLogic

www.precisionbiologic.com  
1.800.267.2796

Frozen diagnostic products. One simple idea, many cool benefits.