

Interlab QC: Chromogenic Factor VIII


Logging in

1. Go to www.precisionbiologic.com and click  (the login icon) located at the top and bottom of each web page or go to my.precisionbiologic.com.
2. Enter your Username (email address) and Password then click **Login**.

Note: If you've forgotten your password, click the **Forgot your password?** link.

3. Click the Interlab QC icon.

Your screen will now look like this:



Using Interlab QC: Chromogenic Factor VIII

1. Click **Enter Chromogenic Factor VIII Data**, **Edit Data**, or **Generate Reports** as appropriate.
 - **Enter Data:** allows you to choose the month for which you are entering information along with your instrument(s), product information and lot numbers, and then enter your raw data. **Ensure to choose the correct unit for your data (% or IU/mL).**
 - **Edit Data:** allows you to adjust any entries you have made.
 - **Generate Reports:** allows you to see the consolidated information for the month chosen based on the raw data that has been entered. **Note:** If you see any discrepancies, try editing the data and regenerating the report. If this does not resolve the issue, please contact us and we will be happy to help.
2. Complete the fields by following the onscreen instructions. *Note: If you are entering or editing data, be sure to click **Save**.*

Questions?

Contact Customer Care (customercare@precisionbiologic.com, 1.800.267.2796 ext. 1), or Technical Support (techsupport@precisionbiologic.com, 1.800.267.2796 ext. 2).