

## CRYOcheck™ dRVVT Reagents



CRYOcheck LA Check™ (dRVVT Screening Reagent) and CRYOcheck LA Sure™ (dRVVT Confirmatory Reagent) are frozen dilute Russell's Viper Venom Time (dRVVT) reagents. These products can be integrated into a lupus anticoagulant (LA) test panel to provide users with enhanced convenience and increased confidence in testing results.

### FEATURES:

- No reconstitution errors
- Easily adapted to manual, semi-automated, and automated instrumentation (protocols for most instrumentation available upon request)
- Unused reagent may be refrozen once
- Packaging and vials are compact and color-coded for easy storage and location in freezers
- Convenient frozen format – ready to use within minutes of thawing
- 48-hour stability once thawed when maintained at 2 to 8 °C
- Expiry dating of 3 years from the date of manufacture when stored at -40 to -80 °C
- Screening and confirmatory reagents are sold separately to eliminate excess stock of either product in your lab

CRYOcheck LA Check and CRYOcheck LA Sure can be used in conjunction with CRYOcheck Lupus Positive Control, CRYOcheck Weak Lupus Positive Control, CRYOcheck Lupus Negative Control, CRYOcheck Pooled Normal Plasma, and CRYOcheck Platelet Lysate.

CRYOcheck LA Check is a sensitive reagent that detects weak LA samples. Reagent contains Russell's viper venom, phospholipids and anti-heparin agents. CRYOcheck LA Sure has a high phospholipid concentration that neutralizes LA and therefore demonstrates the correction needed to report positive LA results with confidence. The result is top quality products which can be used as part of a complete LA workup.

	DESCRIPTION	CATALOG #	FORMAT
	LA Check	CHK-10	25 x 1.0 mL
	LA Sure	SUR-10	25 x 1.0 mL



**PrecisionBioLogic**

www.precisionbiologic.com  
1.800.267.2796

Frozen diagnostic products. One simple idea, many cool benefits.